Nikon A1R. I Widefield Microscopy



Nikon A1R. I Widefield Microscopy

System	Nikon A1R confocal
Location	KBC Building 6 th floor Department of Medical Biochemistry and Biophysics
Room	A6 44 15
contact	irene.martinez@umu.se
Price	50kr/hr. Driver License 1000kr

Use this phrase to mention the microscope in your publications: Nikon Eclipse Ti-E inverted microscope with an DU897 ANDOR EMCCD camera controlled by Nikon NIS Elements interface, equipped with:

Objectives

Nikon CFI Plan Apochromat 10x (N.A 0.45) DIC objective

Nikon CFI Plan Apochromat 20x (N.A 0.75) DIC objective

Nikon CFI Plan Apochromat 60x oil (N.A 1.40) DIC objective

Nikon CFI Plan Apochromat TIRF 100x oil (N.A 1.49) DIC objective

Filter Sets in the turret:

Empty position to perform confocal microscopy

DAPI. Excitation: 340-380 nm Emission: 435-485nm

FITC. Excitation: 465-495 nm Emission: 515-555nm

Cy3. Excitation: 512-552 nm Emission: 565-615nm

TxRed. Excitation: 540-580 nm Emission: 600-660 nm

Camera (Quad Tirf filter)

Switch on the system



1- Turn on Switch #1



2- Turn on Switch #2



If you want to have temperature and CO2 follow step #2 if you don't need the incubation system go to step #12

Switch on the temperature and CO2 at least 1h before starting acquisition:

3- Switch on switch #3







4- Turn on CO2 key located in the wall close to the room door.









5-. Make sure the objectives are placed all the way down. ALWAYS lower the objective to avoid any lense damage. Place the stage incubator carefully





6- Screw the chamber. Make sure you are in the right positions

7- Add 3ml MiliQ water to the reference dish.

8- Put the chamber lid









9- Switch on the temperature from the back

10- Switch on the Air Pump





11- Make sure the mix is correct and stable:

For 5%CO2:

0.8 l/min Air 0.04 l/min of CO2

Or

0.6 l/min Air 0.03 l/min of CO2



If not, please contact BICU personnel: irene.martinez@umu.se



12- Switch on the Epi-Fluorescence Lamp

13- Wait few minutes and then: Switch on the Computer

14- Double click Nis-Elements AR 4.20.01









15- Login with your Name from the list and type your Password

Login X User name: Password: Login Close

16- Select ANDOR camera

NI	S-Elements AR 4.20.01 (Build 982) 64bit - Driver selection	
	Nikon Confocal	•
	No Grabber	
L	+ ANDOR	



Image Acquisition



Main menu:

🕞 😪 🦲 🔞 🔞 🖉 👻 🖉 SH5 💿 AOTE 🗛 EPI 🛀 1.00x 👻 🖑 🦣 Cy3 Eyepiece DAPI 📶 Texas Red DIC TIRE

1- Right click and select the following panels: Ti Pad, DU897 Settings, ND acquisition







2- Select the objective with the TiPad menu or with the microscope botton on the left side

- 3- Add oil to the objective if needed (60X oil 100X oil)
- 4- Place the specimen:

1.00x

BICU

5- Press Eyepiece from the Optical Configuration menu

NME

The last settings will be loaded automatically (blue arrows). In order to change filters (to visualize different wavelengths), use again the microscope or the TiPad menu. Press shutter if is not open

Cy3 Eyepiece DAPI FITC Texas Red DIC

TIRF







6- Focus the specimen.

movement accuracy adjustment

Move the stage with the joystick







7- Press the OC

8- Press Live from the main menu



Cy3 Eyepiece DAPI FITC Texas Red DIC

9- Adjust Camera Settings: Readout Mode: Readout Speed from 1MHz up to 10MHz EM Gain Multiplier: Conversion Gain:

Bining: Binning allows charges from adjacent pixels to be combined. You will get faster readout speeds and improved signal to noise ratios but reduces spatial resolution.

10- Make a picture: Press camera menu. Save the image as nd2 file.



from the main

TIRF







ND Acquisition: With this menu you will automatically acquire images and the Options listed bellow

Pathway to Save your files _

Options

Time: Time lapse experiments.XY: Image different positions.Z: Z Stack

 λ: Optical Configurations to perform an automatic acquisition with different wavelenghts.
 Large Image: Acquisition and stiching

of a large image.





How to build an experiment with different OC. Select the OC from the main menu:



Adjust camera settings for that OC . An exclamation mark will show up. To keep those settings: right click and select: "Assign current camera settings"





Go to ND Acquisition menu and select \square λ



Choose the Optical Configuration you want to use from the list and load them within the SetUp.

To take an image you just have to press Run Now and the camera will acquire the OC selected and the files will be under the pathway you have previously selected.

In order to perform a complex experiment you can select the other parameters within the Set Up and add diferent positions to acquire.



ND Acquisition ×				
Experiment: ND Acquisition				
λ:				
Save to File				
Path: G:\A1 Image:	s Continued\Irene\L	.arge image	Browse	
Filename: conv 4x4.nd2	2		Record Data	
Order of Experiment 👻 T	iming			
🔲 🏶 Time 🔲 🏭 XY 🛛] ∰ z 🗹 🗳 λ	🔲 루 Large Image		
Setup		🕈 Add 🚮 🗐	XX	
Optical Conf.	Name	Comp. Color	Focus Of	
Evepiece	Eyepiece		-	
<no configuration=""> Eyepiece DAPI FITC Cy3 Texas Red DIC TIRF <define new=""></define></no>				
Close active Shutter d	uring Filter Change			
Use Rabo Define Ra				
			Advanced >>	
Load 👻 Save 🔻 🛛	lemove* 🥳	Run Z Corr 1 time lot	op 🥳 Run now	

Time Lapse:

Go to ND Acquisition menu and select Time

You can choose the Interval between acquisitions and the duration. The software will calculate the Loop number

Press Run Now to start the experiment





Positions:

Go to ND Acquisition menu and select
Positions

Press Add Add to add a new position. The coordinates will be kept, then you can move to another area an select a new position. You can travel within the positions selected by clicking with the mouse on each position (blue arrow). Make sure "Move Stage to selected Point" is selected $\Rightarrow \square$

If you are using oil objectives do not travel much along the xy to avoid lensdamages.





Large Image:

Go to ND Acquisition menu and select ☑

Large Image. To use this option you have to select λ too.







Z Stack: Go to ND Acquisition menu and select \square Z.

Pathway to save and name the file

Zstack Top position.Objective is close to the slide

Zstack Botton position. Objective is close to the slide

Step thickness Bottom and Top position







ND sequence acquisition: You can run a sequence on ND acquisitions to perform complex experiments







ND Sequence Acquisition: Time Lapse, Zstack and Time Lapse

Pathway to save the file

Define the ND Sequence. In this case a Time Lapse. Press ok and move to the next Action

				-X		
ND Sequ	ence Acquisiti	ion ×				
Path:	G:\A1 Images	G:\A1 Images Continued\Irene\ND sequence acquisition Browse				
Prefix:	time lapse zst	ack and time I				
Tir	nelapse					
Sequer	nce Definition —		÷	XX		
A	ction	Descriptio	n 🗸			
#1 N	D Acquisition	* ND Acquis	ition: Time(5 loops)	Define		
		ND Sequence / Experiment: T: Order of Experime Time schedul	Acquisition ND Acquisition riment Timing	Record Data	Mage	
				💠 Add 🚽	×	<u> 8</u>
		Phase	Interval	Duration	r Loops	-
		#1	30 sec	▼ 2 min	▼ 5	
Merg	e ND files if pos					
Load	d ▼ S	Close Act	ve Shutter when Idle ansmitted Illuminator off	when Idle (1.00 s)	ne Measurement (0 Ri Use HW seq ents Advanced	OIs) uencer I >>
		Load 🔻	Save 🔻 Remove*		OK Car	ncel





In this case the second Action is a z-stack. We have to load settings from a previous ND file, press ok and then move to the next Action





The las Action is a Time Lapse. Once you have finished mark "Merge ND Files if possible" and Press Run

			×.	
ND Sequence Acquisition ×				
Path: G:\A1 Images Continued\In Prefix: time lapse zstack and time	rene (ND sequence I	acquisition Brow	se	
Timelapse				
Sequence Definition		÷ .		
Action Descr #1 ND Acquisition * ND Ac #2 ND Acquisition * ND Ac #3 ND Acquisition * ND Ac	iption cquisition: Time(5 k cquisition: ZStack(cquisition: Time(5 k	oops) Defii 10) Defii oops) Defii	ne ne ne	
	ND Sequence A	cquisition		X
Experiment: ND Acquisition T:				
	🔽 😫 Time 📋		🖉 A. 🔲 🗮 Large	Thiege
	Time schedule		n 🖗 Add	2 🔿 🗙 🛪
☐ Merge ND files if possible Load ▼ Save ▼	Phase	Interval 30 sec	Duration Z min	
	Close Activ	ve Shutter when Idle nsmitted Illuminator off w	Perform Ti hen Idle (1.00 s) E	me Measurement (0 ROIs)
	Load 🔻	Save 🔻 Remove*		OK Cancel



On going experiment...

Description		ND Sequence Acquisition × Path: G: \A1 Images Continued\Irene \ND sequence acquisition Browse Prefix: time lapse zstack and time I Browse	ND Progress Experiment overall progress: T:	X
		Timelapse Sequence Definition Action Description #1 ND Acquisition * Time(5 loops) #2 ND Acquisition * ND Acquisition: Time(1 loop), ZStack(3) Define #3 ND Acquisition * ND Acquisition: Time(5 loops)	Time elapsed: 0:00:56 Time remaining: 0:01:07 Experiment Status: • Time (Equid.) 3 of 5, (1: 30.04s, Step: 1121.41s, Fps: 0.00) configuring Remaining disk space : 320GB Detail Info Next loop starts in 3 sec.	
307 3140 4		✓ Merge ND files if possible Load Save Remove Finish	Events Next Loop Next Priase Start Priase: #1 Pause & Refocus Finish Abort	
Form 5 10 15 20 25 30 T > EGFP 0.42 µm/px Mono 12bit: 512 x 512 pixels	35			



This is the final image with the nd files merged





SWITCH OFF THE SYSTEM

1- Lower the objectives and remove the excess oil from the objective with lens paper

Don't switch off the microscope if someone has booked within 1 hour after you.

2- ALWAYS Check if that person is going to use the microscope and let him/her know the microscope will stay on.

If you are the last user or next is comming in more than 1 hour:

3- Close Nis Elements

4- Transfer your files: connect the computer to internet to transfer your folder from: Computer / Iomega HDD 3TB (G:) A1 Images Continued to the server. There is a link to the server in the desktop. Once the transfer is done, please unplug the computer from the network for security reasons.

It is forbidden the use of usb or hard disks to transfer your files. You cannot use the computer, hard disk or the server to keep your data as a storage place The server is only to transfer your files. 5- Shut down the computer

6- Turn off the Lamp and check the running time. Write down that info in the logbook.





If you have used the incubation follow the steps bellow, if not, go to step #11

7- Turn off the Temperature

8- Turn off the Air pump



9- Turn off switch #3

10- Close the CO2 main key in the wall





11- Turn off Switch #2 and #1



- 12- Put the plastic cover to protect the microscope
- 13- Sign the logbook with all the information

If you have noticed any anomaly during the acquisition please contact Irene Martinez:

irene.martinez@umu.se

Cell phone: 072 220 6774 ONLY FOR URGENT MATTERS